# Total Lung Capacity (TLC) MiniBox<sup>™</sup> vs. Body Plethysmography

Study conducted at the VA Medical Center, Cincinnati, OH

## INTRODUCTION

The MiniBox system is based on a novel technique for the measurement of lung volumes. The system is designed to measure lung volumes without the utilization of gases or a Plethysmographic cabin. The system sits on a desk and allows a rapid and user-friendly measurement.

The MiniBox technique has been validated in a prospective multi-center study on over 134 subjects, in Israel.

The purpose of this study, as well as others to come, is to test the accuracy of the measurement on different populations, in different geographies, against several Body Plethysmography devices.

## METHOD


This study was held at the VA Medical Center in Cincinnati, and included 50 subjects. Each of the patients enrolled was measured at the same day in random order on the MiniBox and the reference Body Plethysmography device, manufactured by MGC Diagnostics.

## DISCUSSION


Our results show that TLC\_MiniBox is remarkably accurate compared to TLC\_Pleth across the entire population studied. Among our prospective cohort of 50 subjects, who were healthy or had varying severities of obstructive and restrictive diseases, TLC\_MiniBox correlated well with TLC\_Pleth (adjusted r2 = 0.85).

Furthermore, the study showed a CV (Coefficient of Variation) of 9.1%, lower than the CV shown in our multi-center study (134 subjects) of 12.1%, or the CV shown by O'Donnell et al (CHEST, 2010) between Helium Gas Dilution and Computerized Tomography (CT) versus Body Plethysmography of 18.9% and 15.6%, respectively.

	No. subjects	50
STUDY DATA	Male / Female	28/22
Body Box device-	Age (years)*	71.1
MGC Diagnostics.	Height (cm)*	166.9
	Weight (Kg)*	75.4
	No. Healthy	18
	No. Unhealthy	32
	* Mean	



**Figure 1:** Associated Bland-Altman plots comparing MiniBoxTM TLC (TLC\_MiniBox) to Plethysmographic TLC (TLC\_Pleth) for all subjects. The solid lines represent the mean bias while the dashed lines represent the upper and lower limits (±1.96\*SD). The mean coefficient of variation (CV) is displayed within each graph.



**Figure 2:** Scatter plots of plethysmographic TLC (TLC\_Pleth) vs. MiniBoxTM TLC (TLC\_MiniBox) for all subjects. For subjects that were measured more than once on the device, the TLC is presented as the average value of all measurements. The dashed lines represent the unity line and the dotted lines represent the confidence intervals. The linear regression equation and the adjusted R2 are displayed within each graph.

### COMPARATIVE STUDIES RESULTS VERSUS BODY PLETHYSMOGRAPHY

	MiniBox multi-center study, in Israel	MiniBox Boston study	MiniBox Cincinnati study	O'Donnel et al. (CHEST, 2010) – Helium Gas Dilution	O'Donnel et al. (CHEST, 2010) – CT (Computerized Tomography)
No. subjects	134	27	50	149	149
CV [%]	12.1%	10.9%	9.1%	18.9%	15.6%

Pulm**&**n

PFT Outside The Box

PulmOne USA Tel: 1-844-PulmOne (1-844-785-6663) Email us: usinfo@pulm-one.com Visit us: www.pulm-one.com © Copyright 2017. All rights reserved. MiniBox is a trademark of PulmOne.

# Total Lung Capacity (TLC) MiniBox<sup>™</sup> vs. Body Plethysmography

Study conducted by Dr. David J. Kanarek, Pulmonologist, Massachusetts General Hospital, Boston, MA

## INTRODUCTION

The MiniBox system is based on a novel technique for the measurement of lung volumes. The system is designed to measure lung volumes without the utilization of gases or a Plethysmographic cabin. The system sits on a desk and allows a rapid and user-friendly measurement.

The MiniBox technique has been validated in a prospective multi-center study on over 134 subjects, in Israel.

The purpose of this study, as well as others to come, is to test the accuracy of the measurement on different populations, in different geographies, against several Body Plethysmography devices.


### **METHOD**

This study was held at a private practice in Boston, and included 27 subjects. Each of the patients enrolled was measured at the same day in random order on the MiniBox and the reference Body Plethysmography device, manufactured by Morgan Scientific, Inc.


### DISCUSSION

Our results show that TLC\_MiniBox is remarkably accurate compared to TLC\_Pleth across the entire population studied. Among our prospective cohort of 27 subjects, who were healthy or had varying severities of obstructive and restrictive diseases, TLC\_MiniBox correlated well with TLC\_Pleth (adjusted r2 = 0.85).

Furthermore, the study showed a CV (Coefficient of Variation) of 10.9%, lower than the CV shown in our multi-center study (134 subjects) of 12.1%, or the CV shown by O'Donnell et al (CHEST, 2010) between Helium Gas Dilution and Computerized Tomography (CT) versus Body Plethysmography of 18.9% and 15.6%, respectively.



**Figure 1:** Associated Bland-Altman plots comparing MiniBoxTM TLC (TLC\_MiniBox) to Plethysmographic TLC (TLC\_Pleth) for all subjects. The solid lines represent the mean bias while the dashed lines represent the upper and lower limits (±1.96\*SD). The mean coefficient of variation (CV) is displayed within each graph.



**Figure 2:** Scatter plots of plethysmographic TLC (TLC\_Pleth) vs. MiniBoxTM TLC (TLC\_MiniBox) for all subjects. For subjects that were measured more than once on the device, the TLC is presented as the average value of all measurements. The dashed lines represent the unity line and the dotted lines represent the confidence intervals. The linear regression equation and the adjusted R2 are displayed within each graph.

#### COMPARATIVE STUDIES RESULTS VERSUS BODY PLETHYSMOGRAPHY

	MiniBox multi-center study, in Israel	MiniBox Boston study	MiniBox Cincinnati study	O'Donnel et al. (CHEST, 2010) – Helium Gas Dilution	O'Donnel et al. (CHEST, 2010) – CT (Computerized Tomography)
No. subjects	134	27	50	149	149
CV [%]	12.1%	10.9%	9.1%	18.9%	15.6%

PulmOne USA Tel: 1-844-PulmOne (1-844-785-6663) Email us: usinfo@pulm-one.com Visit us: www.pulm-one.com © Copyright 2017. All rights reserved. MiniBox is a trademark of PulmOne.

